QED MEDICAL



ALTAIR HI 2 FIBER OPTIC HEADLIGHT

Brilliant line of sight intensity, maximum control and the finest engineering and components brought together to provide the optimal surgical lighting experience.

FEATURES:

Comfortable, adjustable headband Clean, consistent illumination Field measures 11.5 cm at a 40 cm working distance Specially coated, three-lens system delivers exceptional focus & intensity

Precision internal iris provides pinpoint spot adjustability Our side-mount fiber optic bundle features flexibile durability to bring comfort, intensity and longevity to your system Powered by QED Medical's Xenon, Halogen or LED Light Source Available with 7' or 10' Fiber Optic Cable in Black or Iridescent Blue

Autoclavable positioning bar provides sterile adjustability Quick disconnect cable offers ease of movement Collar and gown clips keep the fiber optic cable secure and contained



ALTAIR SA FIBER OPTIC HEADLIGHT

Small, light, brilliant: The Altair SA Fiber Optic Headlight delivers exceptional intensity in a sleek unobtrusive system.

FEATURES:

Comfortable, adjustable Standard or Flex Headband
Clean, consistent illumination
Simple, lightweight clip-on system for easy integration with
popular magnification loupes
Precision iris adjustment provides a clean, even, variable spot
Pivot points allow for exceptional coaxial alignment
Component design for low cost, simple maintenance
Powered by QED Medical's Xenon, Halogen or LED Light Source
Available with 7' or 10' Fiber Optic Cable
Excellent intensity from examination to the OR
Quick disconnect cable offers ease of movement
Collar and gown clips keep the fiber optic cable secure
and contained



Nova 300-watt Xenon



Nova 200-watt Halogen



Beamer 10-watt LED

WWW.QED-MEDICAL.COM SALES@QED-MEDICAL.COM 750 ENTERPRISE DRIVE, LEXINGTON, KY 40510 800-513-2256 859-231-0338 ph 859-231-0376 fpx



