



**THE C3 3-WATT LED PORTABLE LED HEADLIGHT SYSTEM**

Bring brilliant, even, adjustable LED intensity into the field with the new C3 LED Headlight from QED Medical.



**FEATURES:**

- 50,000 hour 3-watt LED
- 4,500 Foot Candles
- 6,000<sup>o</sup> k (3,200<sup>o</sup> k LED optional)
- Variable Intensity Control
- Battery Life of 7.5 to 100 hours
- No Battery Memory
- Pinpoint Spot Adjustment
- Lightweight, comfortable, headband options
- Sterilizable Positioning Bar
- 100v - 240v Charger
- Transformer option available

**SYSTEM INCLUDES:**

- C3 Headlight
- One Rechargeable Battery Pack
- Charger
- Carrying Case Optional*



*Optional Headband Mounted Battery Pack for Wireless Connectivity Available*

**CONVERSION KIT:**

Convert the QED Medical, High Q or BFW brand Halogen Headlight Systems easily and affordably to LED with the C3 Conversion Kit:

***Centauri, Centauri 2  
Halogen II, Halogen III***

*Conversion Kit includes LED Light & Harness, Battery Pack & Charger*

*Converting to the LED is as easy as changing a bulb*



*Conversion Kit*

**QED MEDICAL**

SURGICAL LIGHTING PRODUCTS SINCE 1971

**THE NAME BEHIND THE SYSTEMS YOU KNOW**



US MADE SINCE 1971

WWW.QED-MEDICAL.COM  
SALES@QED-MEDICAL.COM  
750 ENTERPRISE DRIVE,  
LEXINGTON, KY 40510  
800-513-2256  
859-231-0338 ph  
859-231-0376 fax