

THE LAWSON LIGHT

The need: deliver high quality intensity, shadow-free, into any body cavity.

The answer: The Lawson Light from QED Medical. When asked by one of the leading Otolaryngologists to integrate the comfort, intensity and versatility of our systems with the specialized needs of the ENT doctor, we responded.

FEATURES:

- QED's original ball and swivel linkage allows the optic to be positioned over the dominant eye
- A tunnel drilled with precision through the optic provides an unimpeded, shadowless, beautifully lit, view deep into the restrictive areas of the body
- Comfortable, adjustable Headband
- Clean, consistent illumination
- Field measures 11.5 cm at a 40 cm working distance
- Specially coated, three-lens system delivers exceptional focus & intensity
- Precision internal iris provides pinpoint spot adjustability
- Our side-mount fiber optic bundle features flexible durability to bring comfort, intensity and longevity to your system
- Powered by QED Medical's Xenon, Halogen or LED Light Source
- Available with 7' or 10' Fiber Optic Cable in Black or Iridescent Blue
- Autoclavable positioning bar provides sterile adjustability
- Quick disconnect cable offers ease of movement
- Collar and gown clips keep the fiber optic cable secure and contained



Available in the configurations:

The Lawson Light is available in the following configurations. Select the System that fits your needs best:

- The Fiber Optic System (Altair HI 2 Headlight)
- Portable LED (The C3 LED System)
- Portable Halogen (The Centauri 2 System)



The C3 LED System



The Centauri 2 Halogen System



Nova 300-watt Xenon



Nova 200-watt Halogen



Beamer 10-watt LED

QED MEDICAL

SURGICAL LIGHTING PRODUCTS SINCE 1971

THE NAME BEHIND THE SYSTEMS YOU KNOW



US MADE SINCE 1971

WWW.QED-MEDICAL.COM
 SALES@QED-MEDICAL.COM
 750 ENTERPRISE DRIVE,
 LEXINGTON, KY 40510
 800-513-2256
 859-231-0338 ph
 859-231-0376 fax