## THE BEAMER

THE BEAMER 10-WATT LED LIGHT SOURCE Provides Illumination for Headlights or Instruments

The Beamer provides beautiful, long lasting, white light to your headlight, retractors and lighted instruments. Lightweight, small, unobtrusive and powerful - the Beamer can sit table top or can be attached to any steel surface with its base mounted earth magnet.

The Beamer from QED Medical: the lean, "green" way to illuminate the field.





Bright:	To 4,350 Foot Candles* / 6000° Kelvin
Long Life:	25,000 Hour Lamp Life
Small:	3.25" x 2.5" x 1.75" Footprint
	82 mm x 64 mm x 44 mm
Light:	Less than 10 Ounces (.28 kg)
Quiet:	Near Silent Cooling Fan
Versatile:	Tri-Port for Maximum Access
Universal:	Operates in 100v - 240v Environments
Mounts to Steel Surfaces Magnetically	
*Altair HI 2 Fiber Optic Headlight @ 12"	





The Beamer 10-Watt LED Light Source powers each of these QED Medical Headlights, the ALOS-TTL Through the Lens Video Headlight System and, via the 3-Port Turret, many other available headlights, lighted instruments & retractors.





SURGICAL LIGHTING PRODUCTS SINCE 1971
THE NAME BEHIND THE SYSTEMS YOU KNOW

WWW.QED-MEDICAL.COM SALES@QED-MEDICAL.COM 750 ENTERPRISE DRIVE, LEXINGTON, KY 40510 800-513-2256 859-231-0338 ph 859-231-0376 fax