

THE CENTAURI 2 PORTABLE HALOGEN HEADLIGHT SYSTEM

Maximize the view during examinations with the versatile Centauri 2 Portable Halogen Headlight System from QED Medical.



2-hour, rechargeable, belt-clip, Battery Pack
110v - 240v Transformer delivers uninterrupted use
110v Charger
(220v/240v Charger not shown)

Comfortable, adjustable Headband
Precision iris adjustment provides a variable spot size
Lens & mirror system delivers a clean, even light into the field
Brilliant 11.5 cm spot at a 40 cm working distance
Autoclavable positioning bar provides sterile adjustability
Pivot points allow for exceptional coaxial alignment
Component design for low cost, simple maintenance
Brilliant 10-watt Halogen Lamp provides 100 hours of use
In-line on/off switch provides easy access control

PORTABLE LAWSON LIGHT

The need: a portable way to deliver high quality intensity, shadow-free, into any body cavity. The answer: The Portable Lawson Light from QED Medical.

When asked by one of the leading Otolaryngologists to integrate the comfort, intensity and versatility of our system with the specialized needs of his field, we responded.



Swivel technology allows the Lawson Light optic to be precisely positioned over the dominant eye for true monocular viewing and exceptional coaxial alignment.

All other features of the Centauri 2 apply to the Portable Lawson Light.



US MADE SINCE 1971

QED MEDICAL

SURGICAL LIGHTING PRODUCTS SINCE 1971

THE NAME BEHIND THE SYSTEMS YOU KNOW

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE,
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax