QED MEDICAL



THE XL 10-WATT LED PORTABLE LED HEADLIGHT SYSTEM

A portable LED Headlight built for the OR: Introducing the XL 10-watt LED from QED Medical.

We've integrated the brilliance of a 10-watt, 6000^o Kelvin LED with our Altair HI 2 optic to deliver full functionality, versatility and intensity into the field in a lightweight and comfortable



system suitable for the OR.

25,000 hour 10-watt LED
6,000 Kelvin
Variable Intensity Control
Battery Life of 4.5 hours on High
Charge Time: 5 hours
LED Battery Life Meter
Easy Access On/Off & Intensity Controls
No Battery Memory
Pinpoint Spot Adjustment

64 mm - 13 mm Spot @ 305 mm 3.5" - .75" Spot @ 16"

2.5" - .5" Spot @ 12"

89 mm - 19 mm Spot @ 406 mm Lightweight, comfortable headband Sterilizable Positioning Bar

100v - 240v Charger

SYSTEM INCLUDES:

XL Headlight Two Rechargeable Battery Packs Charger Carrying Case



XL

Driver Control

Battery Pack



LED Battery Life Meter



QED-7870 XL Headight on Standard Headband

QED-7871 XL Headlight on Standard Headband w/Extended Link

QED-7875 XL Rechargeable Battery Pack

QED-7876 XL Battery Charger

QED-7877 XL Replacement Wire Harness

QED-4452 System Carrying Case

QED MEDICAL

WWW.QED-MEDICAL.COM SALES@QED-MEDICAL.COM 750 EMTERPRISE DRIVE, LEXINGTON, KY 40510 800-513-2256 859-231-0338 ph 859-231-0376 fox

