

# QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

## the XL 10-Watt Portable LED Headlight System

Untethered OR Level Intensity

Variable Intensity Control

Lightweight & Comfortable

Comfortable, Adjustable Headband

Battery Powered

Clean, Even field of Intensity

Lithium Ion Rechargeable Pack

Precision Spot Adjustment

2 Battery Packs Included

4.5 hours Each on Maximum Intensity

Built-in Belt Clip

Removeable, Sterilizable Positioning Bar

LED Readout Displays Battery Charge Cycle

Audible Tone Alert for Battery Change

5-hour Charge Time

Full Motion Optic

### *The XL System includes:*

QED-7870 XL 10-watt Portable LED Headlight  
on Standard Headband

QED-7875 XL Rechargeable Lithium Ion Battery Packs

QED-7876 XL Battery Charger

QED-4452 System Carrying Case

Comfortable

Lightweight

Brilliant

Mobile

OR Ready



*Battery Pack*



*Driver*



*System Carrying/Storage Case*

Shown with  
optional  
Positioning  
Bar

WWW.QED-MEDICAL.COM  
SALES@QED-MEDICAL.COM  
750 ENTERPRISE DRIVE  
LEXINGTON, KY 40510  
800-513-2256  
859-231-0338 ph  
859-231-0376 fax