

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

Exceptional Intensity in a Portable LED Headlight

Ideal for ENT, Oral & Maxillofacial Surgery,
Plastic & Reconstructive Surgery and other Medical
and Dental applications

Untethered, Lightweight & Comfortable

Available with 6000k White Light or 3200k Soft
Light Output

Battery Powered

- Lithium Ion Rechargeable Pack
- 7.5 hours on Maximum Intensity
- More than 100 hours at Low Intensity
- Built-in Belt Clip
- Intensity Control Knob
- 3-hour Charge Time

Variable Intensity Control

Comfortable, Adjustable Headband

Clean, Even field of Intensity

Precision Spot Adjustment

Removeable, Sterilizable Positioning Bar

Full Motion Optic

The C3 System includes:

- QED-9700 C3, 3-watt Portable LED Headlight
on Standard Headband
- QED-9745 C3 Rechargeable Lithium Ion Battery Pack
- QED-9746 C3 Battery Charger

C3 Conversion Kit

Upgrade your QED or BFW Portable
Halogen Headlight System to new LED
Technology with this simple, affordable
Conversion Kit including: LED Wiring
Harness, Lithium Battery Pack and Charger

- QED-9701 C3 Conversion Kit
LED Wire Harness, Lithium Battery
and Battery Charger included



US MADE SINCE 1971

the C3 Portable LED Headlight System

Comfortable

Lightweight

Brilliant

Mobile

Adjustable



WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax