

the EOS
Portable LED Headlight

Comfortable

Lightweight

Brilliant

Mobile

Affordable

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

Extremely Small, Lightweight & Unobtrusive

Exceptional Intensity in a Small Portable Unit

Available on Standard or Flex Headbands

Clips Available for Popular Loupes, Safety & Eyeglasses

Battery Powered

Lithium Ion Rechargeable Pack

7.5 hours on Maximum Intensity

More than 100 hours at Low Intensity

Built-in Belt Clip

Intensity Control Knob

3-hour Charge Time

Excellent for Medical, Dental & Veterinary Uses

Composite Sensitive Filter available for Dental Use

Full Motion Optic

The EOS:

QED-4600 on Standard Clip

QED-4604 on Standard Headband

QED-4609 on Flex Headband

QED-4610 Rechargeable Lithium Ion Battery Pack

QED-4611 Battery Charger

*Also available to Accommodate most Popular Loupes,
Safety & Eyeglasses*



Battery Pack



Flex Headband

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax