

QED MEDICAL

US MADE SURGICAL LIGHTING SINCE 1971

the XL 10-Watt Portable LED Headlight System

Untethered OR Level Intensity

Variable Intensity Control

Lightweight & Comfortable

Comfortable, Adjustable Headband

Battery Powered

Clean, Even field of Intensity

Lithium Ion Rechargeable Pack

Precision Spot Adjustment

2 Battery Packs Included

4.5 hours Each on Maximum Intensity

Removeable, Sterilizable Positioning Bar

Built-in Belt Clip

LED Readout Displays Battery Charge Cycle

Full Motion Optic

Audible Tone Alert for Battery Change

5-hour Charge Time

The XL System includes:

QED-7870 XL 10-watt Portable LED Headlight
on Standard Headband

QED-7875 XL Rechargeable Lithium Ion Battery Packs

QED-7876 XL Battery Charger

QED-4452 System Carrying Case

Comfortable

Lightweight

Brilliant

Mobile

OR Ready



Battery Pack



Driver



System Carrying/Storage Case

WWW.QED-MEDICAL.COM
SALES@QED-MEDICAL.COM
750 ENTERPRISE DRIVE
LEXINGTON, KY 40510
800-513-2256
859-231-0338 ph
859-231-0376 fax